

Diagnostic Equipment

Aneroid Manometer



Speculum Holder



Sure Temp



Ophthalmoscopes



Innovative solutions that streamline workflow, improve patient safety and enhance patient care. Amico Accessories provides mounting solutions for all diagnostic equipment, assuring secure mounting and ergonomic solutions everywhere in the healthcare facility.



Features:

- More flexible positioning over the traditional wall board
- Equipment can be locked down or easily moved
- Compatible with all Welch Allyn 767 series
- Saves time on repairs, replacements and installation
- Can be mounted to any headwall, power column, flatwall, equipment boom or direct rail option

Mounting Options

- A** Amico Horizontal Rail Mount
(Compatible with Amico and Modular)
- V** Amico Vertical Rail Mount
(Compatible with Amico and Modular)
- F** Fairfield Rail Mount (push button optional)
- E** DIN 25mm Rail Mount (30mm, 35mm optional)
- C** 1 1/4" Pole Mount
- M** Monitor Rail Mount
- D** Dovetail Rail Mount
(Compatible with Ohmeda Rail)
- G** Hill-Rom Vertical Rail Mount
- S** HSI Mount-Vertical Rail Mount
- W** Direct Wall Mount

Amico Accessories includes all the adapters required for the above rail systems

Equipment **Mounts**

Amico BP **Aneroid**



- | | |
|--------------------|------------------|
| A AM-CB-A | G AM-CB-G |
| V AM-CB-V | S AM-CB-S |
| F AM-CB-F | |
| E AM-CB-E25 | |

System includes; gauge, inflation system, adult cuff, and cuff basket


Cuff **Basket**



- | | |
|-------------------------|--------------------------|
| A AM-CUFBASK-A | C AM-CUFBASK-R125 |
| V AM-CUFBASK-V | M AM-CUFFBASK-MON |
| F AM-CUFBASK-F | G AM-CUFBASK-G |
| E AM-CUFBASK-E25 | S AM-CUFBASK-S |

Mounting plate prevents cuffs from easily falling out


Welch Allyn 767 **Aneroid**



- | | |
|------------------------|-------------------------|
| A AM-767MNT-A | C AM-767MNT-R125 |
| V AM-767MNT-V | M AM-767MNT-MON |
| F AM-767MNT-F | G AM-767MNT-G |
| E AM-767MNT-E25 | S AM-767MNT-S |

Mounting plate only, BP system not included


Welch Allyn 767 **Ophthalmoscopes**



- | | |
|-------------------------|--------------------------|
| A OPHTHALMNT-A | C OPHTHALMNT-R125 |
| V OPHTHALMNT-V | M OPHTHALMNT-MON |
| F OPHTHALMNT-F | G OPHTHALMNT-G |
| E OPHTHALMNT-E25 | S OPHTHALMNT-S |

Mounting plate only, Transformer not included

Welch Allyn 5091 **Aneroid**



- | | |
|-------------------------|--------------------------|
| A AM-5091MNT-A | C AM-5091MNT-R125 |
| V AM-5091MNT-V | M AM-5091MNT-MON |
| F AM-5091MNT-F | G AM-5091MNT-G |
| E AM-5091MNT-E25 | S AM-5091MNT-S |

Mounting plate only, BP system not included


Welch Allyn Sure Temp or KleenSpec **Holder**



- | | |
|----------------------|-----------------------|
| A SPECHLD-A | C SPECHLD-R125 |
| V SPECHLD-V | M SPECHLD-MON |
| F SPECHLD-F | G SPECHLD-G |
| E SPECHLD-E25 | S SPECHLD-S |

Mounting plate only, devices not included


Trimline 4300 **Aneroid**



- | | |
|-------------------------|--------------------------|
| A AM-4300MNT-A | C AM-4300MNT-R125 |
| V AM-4300MNT-V | M AM-4300MNT-MON |
| F AM-4300MNT-F | G AM-4300MNT-G |
| E AM-4300MNT-E25 | S AM-4300MNT-S |

Mounting plate only, BP system not included

WA Baum Wall or ADC 750W **Aneroid**




- | | |
|-------------------------|--------------------------|
| A AM-750WMNT-A | C AM-750WMNT-R125 |
| V AM-750WMNT-V | M AM-750WMNT-MON |
| F AM-750WMNT-F | G AM-750WMNT-G |
| E AM-750WMNT-E25 | S AM-750WMNT-S |

Mounting plate only, BP system not included

Please visit www.amico.com for all technical specifications, installation instructions, detailed AutoCAD drawings and more Instrument Mounting Solutions.



Distributed By:



SP-DE091 PRINTED JULY 2009

the right **connection**

Toll Free: 1.877.264.2697 • Tel: 905.763.7778 • Fax: 905.763.8587 • Website: www.amico.com